

Name: First Aid Cabinet

Code: 81145

First Aid Equipment

Description

Lockable carbon steel cabinet. White. Equipped with increased first aid supplies required by D.M. 388 15/7/03 attachment 1 and D.L. 81 9/4/08 art. 45. Suitable for working environments with three employees or more. Dimensions: 40x58x20 cm.

Content:

- 1 copy D.M. 388 15/7/03
- 3 packs cotton wool
- 1 bottle antiseptic solution, 250 ml
- 1 bottle hydrogen peroxide, 100 ml
- 1 Plastosan box of 100 assorted plasters
- 1 Plastosan box of 100 plasters 7x2 cm
- 3 haemostatic lacets
- 1 pair of scissors 14.5 cm DIN 58279
- 2 adhesive plaster rolls 5 m x 2.5 cm
- 1 PIC 3 case with: 3 liquid soap bags, 3 antiseptic wipes,
- 2 ammonia wipes
- 10 packs of 25 sterile compresses 10x10 cm
- 6 packs sterile compress 18x40 cm
- 2 triangular TNT blankets 96x96x136 cm
- 1 elastic bandage 7 cm
- 2 sterile blankets 40x60 cm DIN 13152-BR
- 2 packs of 10 tissues
- 2 instant Ice Packs, single use
- 1 isothermal blanket, gold/silver 160x210 cm
- 1 pack of 8 assorted bandages
- 1 Emocontrol compression bandage
- 1 tubular net gauze
- 5 sanitary waste bags
- 1 mask with face shield
- 1 clinical thermometer with case
- 2 sterile forceps
- 1 complete kit Pinocchio+Vento
- 3 sterile physiological solution, 500 ml
- 2 bottles antiseptic solution, 500 ml, 10% Povidone-iodine
- 1 Personal sphygmomanometer with manual sthethoscope
- 5 pairs sterile gloves
- 1 First aid multilanguage instructions.



Category	Colours White	Available Sizes
Quantity/Carton 1	Important Information Read carefully the user sheet included. Make sure the item is in good condition, undamaged before use. Once you become aware of a damage, you should immediately replace it.	Recommended for Suitable for working environments with three or more employees. For big communities, companies and industrial plants, restaurants, schools, hotels and holiday villages.
Dexterity	EC Standards	Storage and cleaning Keep in a dry place, far from heat sources.